



张家港麦迪卡医疗器械有限公司  
ZHANGJIAGANG MEDI MEDICAL EQUIPMENT CO., LTD.



CE ISO9001 ISO13485



**Model: YA-3008C HYDRAULIC OPERATION TABLE**

**Size:**(2050±50mm)\*(480±20mm) \*(750-950)±50mm

**Application:** For various surgical operations(head, neck etc.) also orthopedic operations with traction frame

**Functions:**

- 1.HI-LOW:(750-950)±50mm
- 2.Head board fold upward≥30°and fold downward≥90°
- 3.Back board fold upward≥80°and fold downward≥10°
- 4.Leg board fold downward≥90° and fold outward≥90°
- 5.Trendelenburg≥25°and Reverse-Trendelenburg≥30°
- 6.Left and right tilting angle≥20°
- 7.Waist board lifting≥120mm
- 8.Legs can be separated for X-ray film, detachable

**Structure:** Boasting double beds, control on both sides to achieve different positions

**Accessories:** 1 pair armrests, 1 pair waist support frame, 1 pair shoulder supporter, 1 pair leg holders, 1pc anesthesia screen, 1 pair footrests, Imported Y style hermetic ring, Hydraulic pump

**Loading capacity:**180kgs

**Packing:** 133\*78\*105cm=1.1M3 G.W.:300kgs

## Product details

